

Cambridgeshire & Peterborough Community short-term loan equipment



(National Geographic, 2023)

Enable! Explore! Enjoy!

All equipment within this booklet is available for short-term loan (up to a period of 3 weeks, unless in exceptional cases) to children and young people living in Peterborough and Cambridgeshire. The children and young people will need to be known to local NHS health services so that we can ensure the equipment is suitable to meet their needs.

PLEASE NOTE: All equipment listed in this booklet has been kindly funded as gift by Yours Clothing Ltd. This means that the loan of this equipment is not a statutory provision, and therefore not eligible for any additional government funding.

To ensure the longevity of the equipment and for families to enjoy for many years to come, at the point of the loan request you will understand and will be agreeing to the terms of our expectations as outlined below.

General expectations for all equipment loans.

- Give a minimum of 5 working days' notice for your request to be fulfilled.
- Loan of equipment becomes the responsibility of the family in receipt during the loan period.
- You will ensure all parts are present upon delivery.
- You will be liable for any damage or loss of equipment, and it is advised that appropriate insurance cover is in place.
- You are expected to ensure that the equipment is kept in a clean and safe condition during the loan period.
- Ensure equipment is only used the purpose it was intended, and as per the manufacturer's instructions.
- You will not delay in returning equipment as per the loan period.

Taking equipment abroad.

- Equipment can be taken abroad, but this must be specified on the request form.
- You must take appropriate measures to ensure equipment is protected in preparation for transportation.
- It is advised that you have appropriate insurance to cover any potential loss or damage to equipment.

The use of third-party couriers to transport equipment.

- You must take appropriate measures to ensure equipment is protected in preparation for transportation.
- You must ensure that they have appropriate insurance to cover any potential loss or damage to equipment.
- Take appropriate measures to ensure the quality of third-party courier's services and their T&C's.

Use of equipment in care settings.

- Settings must identify a responsible person and provide that person's name and contact number for communication purposes.
- You must ensure that they have appropriate insurance to cover any potential loss or damage to equipment.

To request loan of the equipment please speak with your therapist AND complete the form below - send via email to;

For children known to therapy services (CPFT) in Peterborough - cpm-tr.childrenotcpft@nhs.net

For children known to therapy services (PCC) in Peterborough - – adminhsdm@peterborough.gov.uk

For children known to therapy services (CCS) in Cambridgeshire - CCS-TR.therapyreferrals@nhs.net

This equipment is appropriately maintained. You are responsible to provide your own risk assessment.				
Date of request:				
Item(s) requested:				
When do you want to loan period to start?				
From:				
To:				
Child's NRS ID number:				
Child's name	Child's date of birth:	Child's address:	Parental responsibility:	Contact Telephone Number
Is the item going abroad? Y/N				
Is the item going to a setting? Y/N				
If going into a setting, who will be the responsible person – NAME:			CONTACT NUMBER:	
<p>Please note: If any of the below has not been completed, your order will be rejected.</p> <ul style="list-style-type: none"> • Have you made the family aware of the expectations outlined above? Y/N • Have you provided a copy of the expectations to the family? Y/N • Does the family understand that by requesting the equipment, they are agreeing to the terms of expectations? Y/N • Does the family understand that any equipment on loan to them is their responsibility, and will be liable for any loss or damage? Y/N 				

Delivery and Collection: This will be through NRS.

Contents

Page	Description
5	HIPPOCAMPE ALL-TERRAIN WHEELCHAIR
6	Molift Smart 150 hoist
7	Aerolight Lifestyle Multi-Folding Lightweight & Portable Scooter Ramp
8	Safe and Sound Safety Sleeper
9	Splashy portable bath seat
10	HTS (Hygiene Toileting System)
11	Corner Seat
12	ESCAPE lite – Manual Wheelchair – Standard
13	Mangar ELK lifting cushion
14	Electric Portable Changing Bench (small and large sizes)
15	Go To Seat (size 1 and size 2)
16	Pippin Traveler Bed

Hippocampe All-Terrain Wheelchair

FREEDOM TO GO – ANYWHERE!!



Robust and extremely stable, this smart sports vehicle is light, easy to assemble and packs down in moments.

Designed by disabled people for disabled people, Hippocampe has been stringently tested – conformity, resistance and endurance – to ensure full compliance with CE safety standards.

*it comes with a padded waist belt, harness (small/medium sizes), a headrest and arm rests.

“Thanks to the hippocampe buggy, I can access things with my family like the water, the beach, and even hiking! My summer holidays can be as fun-filled as I want!”



Frame Size	SMALL
Seat Depth	38cm
Backrest height	36-42cm
Seat Width	42cm
Max Occupant Weight	120kg
Weight of buggy	17kg

Useful Information: <https://www.delichon.uk/hippocampe>



Molift Smart 150 Light, foldable, and transportable without tools



Compact and easy to transport

The design of the Molift Smart 150 makes it very easy and flexible to transport when travelling. The compact size also makes it easy to store at home.

Foldable without tools

Molift Smart 150 is designed to be easy and convenient to use at home and to transport when travelling. Being easy to fold without tools is a big advantage.

[f2bbb693/#videos-anchor](https://www.etac.com/en-gb/uk/products/patient-handling/mobile-hoists/molift-smart-150-f2bbb693/#videos-anchor)

“I need assistance with transfers from the floor, so camping has always been very tricky. The Molift hoist enabled me to go camping with my family as it made all my transfers ten times easier!”

Useful Measurements

Item number	M10200
Battery	NiMH 14,4 V-2,6 Ah
Battery charge time (hours)	3
Hoisting height (mm)	1665
Hoisting height min (mm)	210
Weight (KG)	26
Width, chassis, legs open inner	925
Width chassis, outer (mm)	660

Useful information :

<https://www.etac.com/globalassets/inriver-never-delete-this/resources/document/usermanual/bm09299.pdf>

Molift videos: <https://www.etac.com/en-gb/uk/products/patient-handling/mobile-hoists/molift-smart-150-f2bbb693/#videos-anchor>

Aerolight Lifestyle Multi-Folding Lightweight & Portable Scooter Ramp

A perfect vehicle and scooter access ramp, the Aerolight Lifestyle Multi-Folding Scooter Ramp is lightweight and portable yet with an excellent weight capacity.

With a unique design, this ramp has two halves which each fold to half their length, making transportation or moving it far easier. This makes it handy to store the ramp when not in use, or to travel with in a vehicle. It also features rigid handles for safety and ease of carrying.

This ramp can be used for both loading empty mobility scooters or other devices into vehicles, and for you to access your home whilst using your mobility device.

Length 5 ft(150cm)



“I can now stay with my auntie thanks to my temporary ramp! She has a high step to get into her house and I have never been able to go further than the front garden but ow I can get into her house and have sleepovers with my cousins!”

Useful information:

<https://www.healthcarepro.co.uk/mobility-aids/wheelchair-and-scooter-ramps-steps/aerolight-broadfold-scooter-ramp-210cm-7ft?sku=S31134>

https://ccshop.sirv.com/pim/prod/78/1938_b06efd93-0c35-49b6-ad20-57943fe8a3c7

https://ccshop.sirv.com/pim/prod/59/1938_faaea009-2421-4c99-bc17-db6adc72c486

Safe and Sound Safety Sleeper



My daughter used to climb out of a cot all the time and would have injuries meaning I never slept out of worry. I found the safety sleeper and now I can sleep during the night as I know she has a safe space to sleep where she will not hurt herself. We take it everywhere we go with her that involves staying overnight!

The Safety Sleeper provides a totally secure, safe, and comfortable sleep solution for children and adults at risk to themselves or others. Prevents wandering, falling out of bed, injury during seizures and is entrapment proof. Sleep well in the knowledge your loved one is safe, comfy and secure in their Safety Sleeper.

- 100% Safe & Secure
- Easy to Assemble in Minutes
- Easy to transport, weighs less than 19kg (37 lbs)
- Wheeled Travel Suitcase for Easy Transport
- Strong, Stable, Light, and Comfortable
- Perfect for Travel, Respite & Overnights with Family
- Calms & Relaxes Children
- Machine Washable

Dimensions

Width (at base): 90cm

Length (at base): 190cm

Overall height (from base of mattress to top of enclosure): 107cm

Headroom: 81-86cm (depending on mattress used)

Recommended mattress height: 20.5-25.5cm (air mattress supplied is for occasional use only). For everyday use, it is recommended that you use a standard single size mattress

Useful information: [Standard Model - Murrays Medical](#)

Splashy portable bath seat



Splashy is a portable bath seat for children with disabilities which makes bath-time more fun and a lot easier for grown-ups.

It is compact and sits low in the bath, so uses less space and less water, whilst being lightweight, with detachable parts, making it easy to clean. A recline feature means that even children with reduced head control or sitting ability can use Splashy. Splashy comes with four interchangeable foam supports which can be used as a headrest, lateral, pommel or anti slide support

“Splashy is not just for the bath or shower, it can be used at the beach, in a baby pool or for messy play activities. We have even seen it being used in the snow and at water adventure parks.”

“My child is able to access and enjoy bath times when we go on holiday or when he stays with relatives because of the splashy bath seat”.

Features & Benefits



Lightweight and Portable

It's a portable, lightweight and supportive seat that works for bathing or messy play wherever you go.



3 or 5- point harness

Splashy adapts to the child's postural support needs, with the option of a 5-point harness or a 3-point pelvic support.



26 Recline Options

Splashy's base unit lets you switch easily between 26 recline positions, from upright (106 degrees) to full recline (140 degrees).



Adjustable Bumper Supports

Splashy grows with the child. The harness and bumper supports are easily attached using toggles, meaning Splashy's support can be adjusted when required.

Useful information <https://www.leckey.com/products/firefly-splashy>

HTS (Hygiene Toileting System)

Available in Small, Medium, and Large

Innovative Features:

- **In the tub:** The tub stand attaches to the small and medium Rifton HTS seats and bases for comfortable in-tub bathing and hair washing, while the reduced footprint accommodates narrow and wide tub bottoms alike
- **On the road:** The portability kit gives users greater independence – and their families greater freedom to travel. Lightweight folding frame and carry bag (for small and medium sizes) make toileting possible while travelling
- **High-Low Base:** the optional high-low base for the HTS allows the height to be easily adjusted between 51cm to 78cm, assisting with manual handling and to provide access to perform a range of different functions.
- **Tilt-in-space:** With its gas-assisted springs, the Rifton HTS tilts 5° forward for easier transfer, as well as 25° back for showering and hair washing
- **Unique seat pad design:** The unique design of the optional seat pad opening – open to the rear for larger users – facilitates clear access for cleaning and hygiene
- **Optional integral skin foam pads:** For the back and seat, provide comfort for users and easy disinfecting for carers. The pads are easily removed and replaced
- **Load bearing footplate:** For easier transfers

“Having this portable toilet frame meant that I could access the toilet at my holiday club.”



Useful information: [Jiraffe Rifton HTS Brochure 2023 1 Compressed](#)

<https://youtu.be/8zhFPPRyGz0>



Corner Seat

The corner seat is available in a Size 2. It has an adjustable back for height and can be used with or without the upper back support. The seat can help to keep your child well positioned, with their hips neatly back into the seat.



	Measurement	Size 2
Age range (approx.)		1 yr – 6 yrs
User weight	Max	35
Total back support height	Min-max	440-620
Lower back support height	Max	235
Hip width	Max	390
Footprint	Width x length	540 x 565

Useful Information

[Jenx Jiraffe Therapy Brochure 2021_1_5DDR](#)

[Corner Seat - Jiraffe](#)

“I have a specialised seating system at home, but mum and dad are not able to bring that with us when we visit granny and grandad. Instead, they bring my corner seat which means we can all still play together while still meeting my needs”.

ESCAPE lite – Manual Wheelchair – Standard

Product Description

Escape Lite Manual Wheelchair

Minimal weight and an easy folding action make the Escape Lite manual wheelchair the perfect choice for any journey involving a vehicle or wherever storage space is a constraint, thanks to the folding frame and backrest. Designed for comfort for both the user and the attendant, this [lightweight wheelchair](#) is a pleasure to use.

Suitable for indoor and outdoor use, this attendant-propelled wheelchair is an excellent product for occasional use. It weighs just 11.5kg and yet affords a good degree of comfort thanks to padded upholstery, detachable footrests, and full-length armrests.

The attendant-propelled brakes are easy to use with a new push-down feature to lock brakes in a parking position. The Escape Lite Manual Wheelchair is available in a range of different colours and sizes from this area of the shop. As a result, you can choose the ideal wheelchair for you as part of our expansive collection of [wheelchairs](#).

“I have had to have surgery on my leg which means I cannot walk at the moment. I know that it will not be too long before I can walk again but I would have missed going out with my friends; having a wheelchair means I can still do this.”



- Front wheels diameter: 200mm (8")
- Rear wheels diameter: 305mm (12.5")
- Seat width: 460mm (18")
- Seat depth: 410mm (16")
- Seat height (front): 515mm (20.25")
- Seat height (rear): 495mm (19.5")

- Width between arms: 420mm (16.5")
- Backrest height above seat: 400mm (15.75")
- Minimum footrest length: 370mm (14.5")
- Maximum footrest length: 530mm (21")
- Armrest height: 210mm (8")
- Wheelchair length: 900mm (35.4")

- Wheelchair width when open: 530mm (21")
- Wheelchair width when folded: 230mm (9")
- Wheelchair height: 920mm (36")
- Wheelchair weight: 11.5kg (25.4lb)
- Maximum user weight 100kg (15.75st)

Useful information [Escape Lite Manual Wheelchair - Buy Online - Complete Care Shop](#)

Mangar ELK lifting cushion



Mangar Elk Emergency Lifting Cushion

The Mangar Elk [emergency lifting cushion](#) is a compact, portable, battery-powered emergency lifting cushion that has been specifically designed to assist a user safely back to their feet after a fall.

The Elk Emergency Lifting Cushion is totally portable and is supplied with its own battery-powered compressor. That means it can literally be used anywhere, indoors or outdoors, making it particularly effective in small rooms or tight spaces where other lifting equipment cannot be used.

User-friendly and designed to be used with the help of a carer, the person that has fallen is gently rolled onto the lifting cushion and helped to a sitting position. This [moving and handling aid](#) is then gradually inflated whilst the carer stands behind them, supporting them until they are in a high enough raised position to stand.

Use of the Elk minimises the risks of accidents or injuries for both the carer and the user, risks that can be particularly high when assisting a heavier user.

The Elk Emergency Lifting Cushion is used and recommended by most UK NHS Ambulance Trusts as a key piece of emergency equipment.

Product Specification

- One Mangar Elk Emergency Lifting Cushion
- Fully portable emergency lifting system
- Easily operated by just one carer
- Reduces the risks of lifting for carers
- Reduces the risk of injury to the user
- Can be used indoors and outdoors
- Ideal for use in small rooms and tight spaces
- Supplied with rechargeable battery powered compressor
- Folds into compact size when not in use for convenient storage
- Cushion height inflated: 565mm (22")
- Cushion width inflated: 570mm (22.5")
- Cushion depth inflated: 500mm (20")
- Deflated size: 620mm x 260mm x 115mm (24.5" x 10" x 4.5")
- Cushion weight: 3.6 kgs (8 lbs)
- Compressor weight: 4.6 kgs (10 lbs)
- Maximum user weight: 450 kgs (70 st)

Useful information: [Mangar Elk Emergency Lifting Cushion - Complete Care Shop](#) <https://youtu.be/A9lt86iayBY>

Fully mobile Changing table

It requires no installation or fixing. It is electrically operated and is height adjustable. Easy to clean padded surfaces.

Small and large sizes available

Measurement	Size 1	Size 4
Code for mild steel powder coated	3121	3124
Code for stainless steel	3141	3144
Table dimensions length (mm)	1400	2000
Table dimensions width (mm)	750	750
Travel height range (mm)	470-950	470-950
Safety guards width (mm)	840	840
Safety guards height (mm)	165	165
Safety guards operational space (mm)	200	200
Castor centres length (mm)	1090	1090
Castor centres width (mm)	700	700
Nominal clearance (mm)	190	190
Max user weight (kg)	227	227
Product weight (kg)	86	100



Useful information: <https://www.smirthwaite.co.uk/product/mobi-changer/>

Firefly GOTo Seat

A safe, supportive, lightweight and portable seat designed to provide additional upper body postural support for children with disabilities aged 1-8 years. Size 1 and Size 2 available. Useful information: <https://www.leckey.com/products/firefly-goto-seat#product-overview>

GoTo Seat	Size 1	Size 2
Base to Top of Headrest	380mm-540mm [14"-21"]	510mm-680mm [20"-27"]
Lateral Support Width	190mm-260mm [7.5"-10"]	210mm-300mm [8"-12"]
Seat Depth	175mm [7"]	244mm [9.5"]
Max. User Weight	15kg [33lbs]	30kg [66lbs]
Age Range	1-3 years	3-8 years



26 Recline Options

Floorsitter accessory has 26 recline positions, from upright [106°] to full recline [140°].



Adjustable Supports

Height adjustable headrest and laterals allow you to reposition your GoTo as your child grows.



Headrest Option

Standard or advanced headrests available.



Lightweight and Portable

Suitable for attachment to most standard chairs in the home, restaurants and classroom by using the adjustable attachment straps. Can also be attached to shopping trolleys.

Pippin Travel Bed

The Pippin Traveller is a safe, portable bed that can be made in different heights to the individual's requirements. It fits around a standard size single mattress or can be made smaller and supplied with an optional mattress.

It comes with its own transport/storage bags and is ideal for use in hospitals, schools, hospices or taking on holiday.

Flame retardant foam padding is covered with an anti-microbial, waterproof fabric.



Dimensions: 1900mm x 900mm 1000mm high sides.

Useful information: <https://www.kinderkey.co.uk/product/pippin-traveller/>